


THE 3001KL-AMX

ANTIMICROBIAL DIGITAL PLATFORM SCALE WITH EXTENDED HANDRAILS

The all-digital 3001KL-AMX supports up to 1000 pounds and provides the optimum quality, reliability, and ease of use. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale features a large non-skid low profile platform, "live" extended handrails that provide extra support for patient stability, and rugged scuff-proof wheels that enable portability. The scale also offers an antimicrobial keypad, handrails, and transport handle to give patients and providers additional safety. The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare, and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

Antimicrobial
Keypad
Handrails
Transport Handle



The scale provides easy integration and seamless communication with today's EMR applications. The 3001KL-AMX includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 3001KL-AMX is your choice for safety and reliability.

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 3001KL-AMX

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent weight lock

Platform Size (w x d x h)

22" x 15 3/4" x 2 5/8"
(559 mm x 400 mm x 64 mm)

Product Footprint (w x d x h)

22" x 34 1/2" x 52"
(559 mm x 876 mm x 1321 mm)

Product Weight

43 lb / 19.5 kg

Portability

2 wheels for ease of movement

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
3001KL-AMX	Antimicrobial digital platform scale with extended handrails, LB & KG
3001KG-AMX	Antimicrobial digital platform scale with extended handrails, KG only
3001KL-AMX-BT	3001KL-AMX with built-in Pelstar® wireless technology
3001KG-AMX-BT	3001KG-AMX with built-in Pelstar® wireless technology
SS-3001KL-AMX	ScaleSurance 2 year extended warranty for 3001KL-AMX series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscapes.com/ConnectivitySolutions



To learn more, call **1-800-815-6615**
or visit homscapes.com/3001KL-AMX

Health o meter®
Professional



* Not available in all countries