



The solid, industrial design of the 3400KG-C includes a stable, stainless steel base and removable, stainless steel weighing pan that can be autoclaved for sterilization.

## Advantages

- ✓ **Stainless Steel Base and Weighing Surface**  
The stainless steel design is easy to clean and ideal for sterile environments.
- ✓ **Rechargeable Battery Power Source**  
Scale is powered by a 6V rechargeable battery with battery charger included with the scale.
- ✓ **High Resolution Measurements**  
Scale provides highly precise and accurate results with refined 1g resolution.
- ✓ **Backlit Display**  
Easy to read from any angle and in any light condition.
- ✓ **Convenient Weighing Pan**  
The 254mm x 250mm / 10" x 9 7/8" stainless steel weighing pan effectively retains fluids and can be removed for sterilization.
- ✓ **Optional Weighing Tray (3400TRAY)**  
197mm x 248mm / 7 3/4" x 9 3/4" tray is useful for weighing items that require a flat surface. Stainless steel tray is removable and can be autoclaved.
- ✓ **Tare Functionality**  
Container weight can be easily deducted to provide accurate net weight of the contents.



(with optional weighing tray)

Optional additional pans and trays are available, for pans order 3400PAN and for trays order 3400TRAY

### Specifications

- Capacity: 5000g
- Resolution: 1g
- Weighing Pan Surface Dimensions (w x d x h):  
254mm x 250mm x 64mm / 10" x 9 7/8" x 2 1/2"
- Optional Weighing Tray Surface Dimensions (w x d x h):  
197mm x 248mm x 16mm 7 3/4" x 9 3/4" x 5/8"
- Scale Base Dimensions (including handles) (w x d x h):  
292mm x 250mm x 133mm / 11 1/2" x 9 7/8" x 5 1/4"
- Product Footprint (w x d x h):  
276mm x 330mm x 203mm / 10 7/8" x 13" x 8"
- Product Weight: 5kg / 11lb
- Display: 16mm (5/8") LCD Display
- Power Source: 100 – 240 V AC battery charger / 6 volt rechargeable battery (included)
- Functions: Zero, Tare, Auto Zero, Auto Off, Low Battery Indicator, Backlit Display
- 2 Year Limited Warranty

### Options & Accessories

Model	Description
3400KG-C	Digital wet diaper / lap sponge / organ scale with pan
SS-3400KG-C	ScaleSurance 2 year extended warranty
3400PAN	Optional 3 pack of stainless steel pans
3400TRAY	Optional 3 pack of stainless steel trays
20-00018	Battery charger (included), item # for replacement or spare