


## THE 553KL-EHR

# DIGITAL PEDIATRIC SCALE WITH EXTRA-WIDE TRAY AND DIGITAL HEIGHT ROD

The 553KL-EHR pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 44 pounds and 20 kilograms, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The easy to read digital height rod makes measuring infant height fast and easy with clear digital measurements, large head and foot paddles, and convenient Hold function. The scale ensures a secure platform, keeps patients calm for accurate readings, and accommodates infants lying down to toddlers sitting up.

With its user-friendly digital interface, the 553KL-EHR scale includes features such as Zero, Tare, Hold/Release, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 553KL-EHR includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 553KL-EHR is your choice for safety, comfort, and reliability.

**120+** LOCAL  
FIELD SPECIALISTS

**100%**  
CONNECTIVITY

**#1** TRUSTED  
BRAND

**100** YEARS  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 553KL-EHR

weigh**easier**®

## TECHNICAL SPECIFICATIONS

### Capacity

44 lb / 20 kg

### Resolution

0.5 oz / 10g

### EMR Connectivity

Via USB, optional Pelstar® wireless technology

### Display

1" / 25 mm, LCD display

### Power Source

6 AA batteries included, 100-240V adapter (ADPT30) not included

### Functions

LB / KG conversion

LB / KG lock out

Zero

Tare

Hold / release

Auto zero

Auto off

Motion detection

### Tray Size (w x d x h)

24 7/8" x 13" x 1 3/8"

(632 mm x 330 mm x 35 mm)

### Product Footprint

24 7/8" x 16 1/2" x 3 7/8"

(632 mm x 419 mm x 98 mm)

### Product Weight

10 lb / 5 kg

### Digital Height Range

13 3/4" - 31 1/2" / 35 cm - 80 cm

Resolution: 1/16" / 0.1 cm

### Height Rod Power Source

2 AAA batteries (included)

### Height Rod Display

LCD, 5/8" / 16 mm digit height

### Portability

**2210CART**- optional rolling cart

### Warranty

2-year limited manufacturer warranty

## ORDERING INFORMATION

MODEL	DESCRIPTION
553KL-EHR	Digital Pediatric Scale with Extra-Wide Tray and Digital Height Rod, LB & KG
553KG-EHR	Digital Pediatric Scale with Extra-Wide Tray and Digital Height Rod, KG only
SS-553KL	ScaleSurance® 2 year extended warranty for 553 series*
2210CART	Rolling scale/ medical device cart for 553 series of scales
ADPT30	Optional adapter for 553 series of scales
PLUGSET30	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer



### 2210CART

Medical Device Cart for 553 Series of Scales

To learn more, call **1-800-815-6615**  
or visit [homscales.com/553KL-EHR](https://homscales.com/553KL-EHR)

**Health o meter®**  
Professional



\*Not available in all countries.