



- Durable aluminum construction
- Installs easily and quickly with simple tools and provided hardware
- Can be mounted on wall or attached to pediatric tray scale, see options chart below

Specifications

- Range of measurement: 7 – 100 cm/3 – 39”
- Division of measurement: 1 mm / 1/16”
- Product Dimensions: 41” / 104 cm (l) x 1.5” / 3.8 cm (w) x 11.25” / 28.6cm (d)

Model	Description
BABYHR	Mechanical baby height rod (includes wall mounting hardware)
522HR	BABYHR that attaches to 522KL & 522KG
553HR	BABYHR that attaches to 553KL & 553KG