

## Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod

ELEVATE-C

ELEVATE

The Health o meter<sup>®</sup> EMRscale

### Includes Connectivity Kit

Elevate includes a connectivity kit that allows the scale to interface directly with Welch Allyn vital signs monitors.



Improve workflow with this top-of-the-line physician office scale. The built-in digital height rod eliminates the need to manually enter height for computing BMI. The 350° swivel head allows for line-of-sight viewing and patient privacy.

## Advantages

- ✓ **EMR Capable**  
Scale interfaces with a computer, EMR software or other electronic device via Health o meter<sup>®</sup> Professional Connectivity Solutions
- ✓ **350° Swivel Display**  
Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy
- ✓ **Built-in Digital Height Rod**  
The easy-to-use, state of the art digital height rod provides measurements from 23 5/8" - 88" / 60 cm - 222 cm
- ✓ **Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost
- ✓ **Calculates Body Mass Index (BMI)**  
Scale functions include the ability to calculate the patient's BMI
- ✓ **Patented Design**  
The ELEVATE's patented design provides a user-friendly and intuitive interface for measuring patient weight, height and BMI with a convenient rotating display head

**Health o meter<sup>®</sup>**  
Professional  
[www.homscales.com](http://www.homscales.com)

See reverse side for additional specifications

# Heavy Duty Eye Level Digital Scale with Integral Digital Height Rod

ELEVATE-C

**Health o meter®**  
P r o f e s s i o n a l



**ELEVATE**  
The Health o meter® EMRscale



US Patent #D667323

## Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Height Rod Resolution: 0.1" / 0.1 cm
- Platform Size (w x d x h):  
14 1/4" x 14 1/4" x 2 5/8" / 362 mm x 362 mm x 67 mm
- Product Footprint (w x h x d):  
14 3/4" x 24 1/8" x 54 3/8" / 375 mm x 613 mm x 1381 mm
- Product Weight: 28 lb / 13 kg
- Digital Height Rod: 23 5/8" - 88" / 60 cm - 222 cm
- Includes connectivity kit to interface with Welch Allyn® vital signs monitors
- Connectivity: USB, optional Pelstar wireless technology
- 1" / 25 mm LCD display, 35° left / right swivel  
(350° with height rod fully retracted) for ease of viewing
- Functions: LB / KG Conversion, LB / KG Lock Out, BMI, Zero, Tare, Hold / Release, Auto Zero, Auto Off
- 2 wheels for ease of movement
- 100 - 240V adapter (ADPT50) included, 6 C-cell batteries (not included)
- 2 year limited warranty

Model	Description
ELEVATE-C	Heavy duty eye level digital scale with integral digital height rod
ELEVATE-KG	Heavy duty eye level digital scale with integral digital height rod, KG only
ELEVATE-BT	ELEVATE-C with built-in Pelstar wireless technology
ELEVATE-KG-BT	ELEVATE-BT scale, kilograms only
SS-ELEVATE-C	ScaleSurance 2 year extended warranty for ELEVATE-C & ELEVATE-KG*
SS-ELEVATE-BT	ScaleSurance 2 year extended warranty for ELEVATE-BT models*
ADPT50	Adapter included, item # for replacement or spare
PLUGSET50	Plug adapters for use outside North America
PELSTARONE-P	Wireless thermal printer

Health o meter® Professional Scales  
9500 West 55th Street USA  
McCook, IL 60525-7110  
Phone: 1-800-815-6615  
Fax: 1-800-654-7330

\* Not available in all countries

ELEVATE-C\_F&B\_20210412

**Health o meter®**  
P r o f e s s i o n a l  
www.homscales.com