


Introducing the Patient Transfer Scale

- A transfer board with a built-in medical scale provides a weight measurement when transferring a patient
- Presents an accurate weight when it is critical for safe medication management for time sensitive patients
- Can be utilized as part of an existing patient transfer process



Weigh critical or immobile
patients instantly and accurately

*Making Healthcare Professionals' Workflow Weigh**asier***[®]

Health o meter[®]
Professional



UK PATENT
Patent number 2543835



US PATENT PENDING



EU PATENT PENDING



CANADA PATENT PENDING





HONG KONG PATENT PENDING

The Patient Transfer Scale: Every Minute Counts



Key Features / Specifications:

- **Large, clear digital weight display**
- **Portable**
- **Wipe clean surface**
- **Capacity:** 550lb / 250kg
- **Resolution:** 1lb / 0.5kg
- **Power:** Internal rechargeable battery
- **Battery life:** 50 hours
- **Weight:** 25lb / 11.4 kg
- **Dimensions:**
71" x 27 1/2" x 1" / 1805mm x 700mm x 30 mm
- **Model Number:**
PTS-1000KL
PTS-1000KL-C

- When managing a patient with a critical condition like a stroke, **Every Minute Counts.** It is imperative to get an accurate weight quickly when the patient is first admitted to the hospital. The Patient Transfer Scale (PTS) does just that.
- All departments can benefit from the Patient Transfer Scale, particularly Stroke Units, the Emergency Room, ICU, and Radiology.
- Learn more about how the PTS works Visit <https://www.homscales.com/products/patient-transfer-scale>
- Get to know the inventor of the PTS Gillian Taylor. Read her success stories.
 **Twitter: @gilliantaym999**
 **Instagram: @gilliantaylorpts**
- The Patient Transfer Scale was developed and invented by an Emergency Care nurse, with the final design based on feedback from healthcare professionals
... For Nurses by Nurses.

Scott Gottman
Director of Business and Product
Development
(708) 297-6582
sgottman@homscales.com

Health o meter® Professional Scales
(800) 815-6615
www.homscales.com

weigh**easier**®