

THE 201HR-3105

MECHANICAL HEIGHT ROD FOR THE 3105 SERIES OF SCALES

The Health o meter® Professional 201HR-3105 is an optional mechanical height rod for the Antimicrobial 3105 series of scales. The durable aluminum height rod withstands frequent use, and rigorous laboratory testing shows that the 201HR-3105 lasts 2 ½ times longer than similar products.

The easy-to-use smooth sliding headpiece and telescoping rod allow the caregiver to read the height measurement quickly and easily. The folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The 201HR-3105 has an extensive measuring range of 2"-90" in 1/16" increments (6 cm-230cm in 1mm increments). The optional mechanical height rod has upgraded brackets making assembly quick and easy.



TECHNICAL SPECIFICATIONS

Measuring Range

2" – 90" / 6cm – 230cm

Graduation: 1/16" / 1mm

Warranty

2-year limited manufacturer warranty

201HR-3105 Compatible Scales

3105KL-AM

3105KG-AM

3105KL-AM-BT

3105KG-AM-BT

ORDERING INFORMATION

MODEL	DESCRIPTION
201HR-3105	Optional mechanical height rod for 3105 Series

Health o meter®
Professional

Health o meter® Professional 201HR-3105

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/201HR-3105**