

THE 3105KL-AM

HEAVY-DUTY ANTIMICROBIAL PLATFORM SCALE

The 3105KL-AM is a robust antimicrobial scale included in a full line of Engineered with a Purpose™ Acute Care products. This heavy-duty platform scale supports up to 1000 pounds (454 kilograms) and provides the optimum in quality, reliability, and ease of use.

The easy-to-use 3105KL-AM improves workflow with its ability to weigh accurately despite the patient's movement. The smooth-rolling wheels and easy-to-grip handle make moving from room to room simple. The scale offers an optional mechanical and market-leading digital height rod. Other options include a kilograms-only scale and built-in Pelstar® wireless technology. This scale is factory calibrated and does not require calibration before use. If you are looking to calibrate your scale, the 3105KL-AM leverages a quick one-step calibration process in LB or KG.

The low-profile stable base and "live" handlebars provide a safe weighing environment, increasing the confidence of both clinician and patient during the weighing process. Other safety features include antimicrobial protection on the head, keypad, transport handle, handrails, body, and kickplate and an innovative design that eliminates places for dirt and contaminants to collect.

Designed to meet the needs of hospitals, acute care facilities, doctor's offices, and clinics. The 3105KL-AM's high-contrast color screen provides a user-friendly interface that includes BMI, hold/ release, reweigh, recall, zero, KG/LB, tare, pretare and more.

The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale. The scale provides easy integration and seamless communication with today's EMR applications. The 3105KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

THERE'S STRENGTH IN OUR NUMBERS.


When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE



Health o meter®
Professional

Health o meter® Professional 3105KL-AM

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

6 D batteries included, 100-240V adapter (ADPT30) not included

Functions

LB/KG conversion
LB/KG lockout
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto-off
Time/date
Variable auto-off time
Audible/mute sound option
Parent/child weighing function
Everlock® - unit measure lock

Platform Size (w x d x h)

19 3/4" x 21 1/4" x 2 1/8"
(502mm x 540mm x 54mm)

Product Footprint (w x d x h)

23 1/4" x 28" x 52 1/4"
(591mm x 711mm x 1327mm)

Product Weight

86 lb / 39 kg

Portability

2 smooth-rolling wheels and easy grip handle make it easy to move from room to room

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
3105KL-AM	Antimicrobial heavy-duty platform scale, LB & KG
3105KG-AM	Antimicrobial heavy-duty platform scale, KG only
3105KL-AM-BT	3105KL-AM with built-in Pelstar® wireless technology
3105KG-AM-BT	3105KG-AM with built-in Pelstar® wireless technology
245EHR-3105	Digital height rod for the 3105 series of scales
201HR-3105	Mechanical height rod for the 3105 series of scales
SS-3105KL-AM	ScaleSurance® 2 year extended warranty for 3105 series*
ADPT30	Optional adapter for 3105 series of scales
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit hom scales.com/ConnectivitySolutions

Optional Digital and Mechanical Height Rods Available



To learn more, call **1-800-815-6615**
or visit hom scales.com/3105KL-AM

Health o meter®
Professional



*Not available in all countries.

3105KL-AM_F&B_20240201