

THE 1110KL

HIGH CAPACITY REMOTE DISPLAY DIGITAL SCALE WITH LARGE BASE

The 1110KL is a remote display digital scale with a high capacity that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The remote display provides versatility with the ability to be mounted on a wall, placed on a table, or handheld. Designed to meet the needs of doctor's offices, hospitals, clinics, and acute and long-term care facilities, the scale's keypad features antimicrobial protection, and the large non-skid low profile platform requires a minimal step-up providing a safe weighing environment.

The scale's high-contrast color screen provides a user-friendly interface that includes BMI, hold/release, reweigh, recall, zero, KG/LB, tare, pre-tare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The 1110KL was designed with a facilities' ever-changing needs in mind. The optional "live" handlebars provide stability and added safety for the patient and caregiver. An optional digital or mechanical height rod can be added to the handlebar accessory to measure height and weight simultaneously.



The scale provides easy integration and seamless communication with today's EMR applications, the 1110KL includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers, it's not surprising that the 1100 series family of scales is the #1 selling platform scale in the US.

20,000+
INSTALLED NATIONWIDE

120+ LOCAL
FIELD SPECIALISTS

100%
CONNECTIVITY

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 1110KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm

Screen Size: 2 3/4" x 1 1/2" High-Contrast

Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion

LB/KG lock out

BMI

Zero

Tare

Pre-tare

Hold/release

Reweigh

Recall

Auto zero

Auto off

Time/date

Variable auto off time

Audible/mute sound option

Parent/child weighing function

Everlock®

Platform Size (w x d x h)

22" x 15 3/4" x 2 5/8"

(559 mm x 400 mm x 66 mm)

Product Footprint - Base (w x d x h)

22" x 23 1/2" x 4 3/4"

(559mm x 597mm x 121 mm)

Product Weight

27 lb / 12 kg

Remote Cable

6 ft / 1829 mm

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
1110KL	High capacity digital platform scale with remote display, LB & KG
1110KG	High capacity digital platform scale with remote display, KG only
1110KL-BT	1110KL with built-in Pelstar® wireless technology
1110KG-BT	1110KG with built-in Pelstar® wireless technology
1110HB	Optional handlebars for the 1110 series, includes 2 wheels for portability
245EHR-1110	Optional digital height rod for 1110 series (must have an 1110HB)
201HR-1110	Optional mechanical height rod for 1110 series (must have an 1110HB)
SS-1110KL	ScaleSurance® 2 year extended warranty for 1110 series*
SS-245EHR	ScaleSurance® 2 year extended warranty for 245EHR-1110*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit hom scales.com/ConnectivitySolutions

To learn more, call **1-800-815-6615**
or visit hom scales.com/1110KL



*Not available in all countries.



1110KL - High capacity digital platform scale with remote display

1110KL with optional handlebars and wheels

1110KL with optional handlebars, wheels and height rod

Health o meter®
Professional