

THE 499HB

HANDLEBAR ACCESSORY FOR THE 499 SERIES OF SCALES

The Health o meter® Professional 499HB is a unique handlebar accessory for the 499KL and is easy to attach to any new scale using the pre-drilled mounting holes and can also be attached to any existing scale using the easy-to-use drill guides. It provides the stability that many patients need without affecting the accuracy of the weight reading.



TECHNICAL SPECIFICATIONS

Product Footprint with 499KL
19 1/4" x 21 1/8" x 45 1/2"
(489 mm x 537 mm x 1156 mm)

Product Weight with 499KL
39 lb / 17.7 kg

Warranty
2-year limited manufacturer warranty

499HB Compatible Scales

499KL
499KG
499KLAD
499KGAD
499KLHR
499KGHR

ORDERING INFORMATION

MODEL	DESCRIPTION
499HB	Optional handlebar accessory for 499 series scales

Health o meter®
Professional

Health o meter® Professional 499HB

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/499HB**